

Data Input

Upload EDF File

chb01_01.edf

40.4MB

+

Upload .seizures File (Optional)

Upload

200MB per file • SEIZURES, TXT

UI Customization

Light Mode

File loaded successfully!

Model Selection

Select ML Model

DA-GRL

Model Requirements: 18 channels

Visualization Parameters

Select EEG Channel to View

FP1-F7

Zoom Window (For UI Visualization Only)

0.00

60.00

Processing Options

Apply Bandpass Filter (0.5 - 45 Hz)

Apply Z-Score Normalization

Pipeline Parameters

Window Size (seconds)

4.00

Step Size (seconds)

2.00

SHAP XAI Options

Enable SHAP Channel Selection

Top-K Channels to Keep

5

XAI SEIZURE PREDICTION APP

Robust EEG seizure detection supporting multiple models with pipeline optimization and SHAP explainability.

Analysis Dashboard

Historical Reports

Time-Series

Frequency

Spectral Analysis

Model Prediction &

Channel Heatmap

Custom Feature

Spectral Topography

SHAP Interpretability

Final Clinical Report

Brain Topography: FFT Relative Power

Spatial distribution of spectral power across the scalp.

Electrode positions synchronized successfully.

Delta (0.5-4Hz)

Theta (4-8Hz)

Alpha (8-13Hz)

Beta (13-30Hz)

Visualization generated using 18 mapped scalp electrodes.

localhost:8501/#running-inference-da-grl

1/1